Description

The HemiCAP® Contoured Articular Prosthetic incorporates an articular resurfacing component and a taper post fixation component that mate together via a taper interlock to provide stable and immobile fixation of the implant and stress bearing contact at the bone/ prosthetic interface.

Materials

Articular Resurfacing Component: Cobalt-Chromium Alloy (Co-Cr-Mo)

Surface Coating: Titanium (CP Ti)

Taper Post: Titanium Alloy (Ti-6Al-4V)

Indications

For the reconstruction of painful and/or severely disabled shoulder joints resulting from post-traumatic degenerative disease or avascular necrosis. The humeral head and neck should be of sufficient bone stock to support loading. The rotator cuff should be intact or reconstructable. The device is a single use implant intended to be used with bone cement.

Patient selection factors to be considered include:

1. Need to obtain pain relief and improve function.
2. Patient age as a potential for early-age-revision of total joint arthroplasty.
3. Patient overall well-being, including ability and willingness to follow instructions and comply with activity restrictions.

Contraindications

Absolute contraindications include:

1. Defects that are not localized.
2. Defects that are located on joint surfaces that are discontinuous.
3. Inflammatory degenerative joint disease, rheumatoid arthritis, infection, sepsis, and osteomyelitis.
4. Patients that have a known sensitivity to Cobalt-Chrome alloys typically used in prosthetic devices.

Relative contraindications include:

1. Uncooperative patient or patient incapable of following preoperative and postoperative instructions.
2. Metabolic disorders which may impair the formation or healing of bone.
3. Infections at remote sites which may spread to the implant site.
4. Rapid joint destruction or bone resorption visible on roentgenogram.

5. Chronic instability or deficient soft tissues and other support structures.

6. Vascular or muscular insufficiency.

Warnings

Improper selection, placement, positioning, alignment, and fixation of the implant components may reduce the service life of the prosthetic components. Inadequate preparation and cleaning of the implant components mating surfaces may result in improper fixation of the device. Improper handling of the implants can produce scratches, nicks or dents that may have adverse clinical effects on mating joint surfaces. Do not modify implants. The surgeon shall be thoroughly familiar with the implants, instruments, and surgical technique prior to performing surgery.

When defining offsets of articular surfaces, care should be taken to ensure that instruments are properly aligned and mated with taper in Taper Post. Visually confirm distal tip of contact probe is making contact on articular surfaces and free from any soft tissue structures to ensure accuracy. Use light pressure on contact probe to slightly indent articular surface at each measurement point, ensuring that the selected implant will be flush or slightly recessed with the articular surface.

When placing implant, carefully trim articular cartilage debris around margin of implant. Remove bone particles and lavage thoroughly. To ensure mechanical interlock of the Taper Post and implant, carefully clean Taper Post taper with provided instruments. All drilling or reaming should be done at slowest speeds possible with vigorous lavage to minimize heat effects to adjacent bone and cartilage tissues.

Accepted practices in post operative care should be used. The patient is to be instructed and monitored to ensure a reasonable degree of compliance to post operative instructions and activity restrictions. Excessive activity, impact, and weight gain have been implicated in the reduction of the benefit and service life of prosthetic devices.

Precautions

HemiCAP[®] implants are intended to be fitted and installed with the HemiCAP[®] instrument set. Use of instruments from other systems may result in improper implant selection, fitting, and placement which could result in implant failure or poor clinical outcome. The HemiCAP[®] instrument set should be regularly inspected for any signs of wear or damage. Do not reuse implants. Reuse of single use devices can increase

the risk of patient infection and can compromise service life and other performance attributes of the device.

Possible Adverse Effects

1. **Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions. Particulate wear debris and mild tissue discoloration from metallic components have been noted in other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.**
2. **Infection or allergic reaction.**
3. **Loosening, migration or loss of fixation of implant.**
4. **Fretting and crevice corrosion can occur at the interface between the implant components.**
5. **Fatigue fracture of the implants as a result of bone resorption around the implant components.**
6. **Wear and damage to the implant articulating surface.**
7. **Wear and damage to the adjacent and opposed articular cartilage surfaces or soft tissue support structures.**
8. **Intraoperative or postoperative bone fracture.**

Sterility

Prosthetic components are sterilized by exposure to gamma irradiation. Do not resterilize. Do not use components if packaging is opened or damaged. Do not use components if beyond expiration date.

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
1. Use **Drill Guide** to locate the axis normal to the articular surface and central to the defect. Choose the correct **Drill Guide** diameter sufficient to circumscribe the defect. Confirm the

appropriate **Articular Component** diameter by matching it to the **Drill Guide** diameter. Place **Guide Pin** into a Cannulated Powered Drill and secure at the etch marking on the **Guide Pin**. Advance **Guide Pin** into bone making sure that it is central to the defect. (It is important to verify the **Drill**

Guide is seated on the curved surface such that four points of contact are established on the articular surface. A normal axis and correct **Articular Component** diameter are necessary for proper implant fit.)

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2. Place **Cannulated Drill** over **Guide Pin** and drive until the proximal shoulder of the **Drill** is flush with the articular surface. Tap hole to etched depth mark on **Tap**. (Use lavage during drilling to prevent possible tissue damage from heat effects.)

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
3. Prior to inserting the **Taper Post**, thoroughly cleanse the pilot hole of any debris and then inject the cement in a retrograde fashion from the end of the hole upwards.

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
4. Place the **Driver** into the **Taper Post** and advance the **Taper Post** until the line on the **Driver** is flush with the cartilage surface.

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
5. Clean taper in **Taper Post** with **Taper Cleaner**. Place **Trial Cap** into **Taper Post** to confirm correct depth of **Taper Post**. The height of the **Trial Cap** must be flush or slightly below the existing articular cartilage surface to avoid the **Articular Component** from being placed proud or above the surface of the defect. Adjust depth if needed using the **Driver** to rotate the **Taper Post** (rotate clockwise to advance and counterclockwise to retract). Remove **Trial Cap**.

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
7. Remove **Centering Shaft** and replace with **Guide Pin**. Advance **Circle Cutter** onto the articular surface by twisting the **Circle Cutter** back and forth avoiding any bending of the **Guide Pin**.

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
8. Choose the appropriate **Surface Reamer** based on the offsets. Confirm selection by matching the color code on the **Articular Component** package with the colored band on the **Surface Reamer** shaft. Drive **Surface Reamer** over **Guide Pin** until it contacts the top surface on **Taper Post**. (Use lavage during drilling to prevent possible tissue damage from heat effects.) Make sure not to bend the **Guide Pin** during drilling as it may result in **Articular Component** malalignment.

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
11. Before placing the **Articular Component** on the **Implant Holder** make sure that sufficient suction is present to hold the device on the distal suction cup. Align the **Articular Component** on the **Implant Holder**. For non-spherical

Articular Components orient the etch marks on the back of the **Articular Component** with the etch mark on the handle of the **Implant Holder**. Align the **Articular Component** with the appropriate offsets. Insert into taper of **Taper Post**.

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
12. Use a slight tap on the **Impactor** to seat **Articular Component**. Progressively tap the **Impactor** until the **Articular Component** is firmly seated on the bone.

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	Consult Instructions for Use
	Caution
	Do Not Reuse
	Sterilized using Radiation
	Sterilized using Vaporized Hydrogen Peroxide
R_x ONLY	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician
	Do Not Use If Package Is Damaged
	Lot Number
	Expiration Date
	Manufacturer
	Catalogue Number
	Single sterile barrier system with protective packaging inside
	Double sterile barrier system
	Date of Manufacture

Arthrosurface, Inc.
319 Manley St.
West Bridgewater, MA 02379
508-520-3003
www.arthrosurface.com

This product is covered by one or more of U.S. Patent Nos. 6,520,964;
6,610,067; 6,679,917; other patents and other patents pending.
HemiCAP[®] is a trademark of Arthrosurface, Inc. U.S.
© 2024 Arthrosurface, Inc. All rights reserved. Printed in the U.S.A.

PN 3001-2000 REV L