  
arthrosurface®

 **HemiCAP**  
Hemiarthroplasty Components  
MTP Instructions for Use

### **Description**

The HemiCAP® Contoured Articular Prosthetic incorporates an articular component and a cancellous taper post component that mate together via a taper interlock to provide stable and immobile fixation of the implant and stress bearing contact at the bone/prosthetic interface.

### **Materials**

Articular Component: Cobalt-Chromium Alloy (Co-Cr-Mo)

Surface Coating: Titanium (CP Ti)

Taper Post Component: Titanium Alloy (Ti-6Al-4V)

### **Indications**

Hemiarthroplasty implant for the metatarsophalangeal joint for use in the treatment of patients with degenerative and post-traumatic arthritis in the metatarsal joint in the presence of good bone stock along with the following clinical conditions: hallux valgus or hallux limitus, hallux rigidus, and an unstable or painful metatarsal/phalangeal (MTP) joint. The device is a single-use implant intended to be used with bone cement.

### **Patient selection factors to be considered include:**

1. Need to obtain pain relief and improve function.
2. Patient age as a potential for early-age revision of total joint arthroplasty.
3. Patient overall well-being, including ability and willingness to follow instructions and comply with activity restrictions.
4. Failure of previous conservative treatment options in correcting deformity and achieving pain relief.

### **Contraindications**

#### **Absolute contraindications include:**

1. Significant bone demineralization or inadequate bone stock.
2. Inadequate skin, musculotendinous or neurovascular system status.
3. Inflammatory or rheumatoid arthritis, infection, sepsis, and osteomyelitis.
4. Patients that have a known sensitivity to Cobalt-Chrome alloys typically used in prosthetic devices.

**Relative contraindications include:**

1. Uncooperative patient or patient incapable of following preoperative and postoperative instructions.
2. Metabolic disorders which may impair the formation or healing of bone.
3. Infections at remote sites which may spread to the implant site.
4. Rapid joint destruction or bone resorption visible on roentgenogram.
5. Chronic instability or deficient soft tissues and other support structures.
6. Vascular or muscular insufficiency.

**Warnings**

Improper selection, placement, positioning, alignment, and fixation of the implant components may reduce the service life of the prosthetic components. Inadequate preparation and cleaning of the implant components mating surfaces may result in improper fixation of the device. Improper handling of the implants can produce scratches, nicks or dents that may have adverse clinical effects on mating joint surfaces. Do not modify implants. The surgeon shall be thoroughly familiar with the implants, instruments, and surgical technique prior to performing surgery.

When mapping articular surfaces, care should be taken to ensure that instruments are properly aligned and mated with taper in Taper Post. Visually confirm distal tip of contact probe is making contact on articular surfaces and free from any soft tissue structures to ensure accuracy. Use light pressure on contact probe to slightly indent articular surface at each measurement point, ensuring that the selected implant will be flush or slightly recessed with the articular surface.

When placing implant, carefully trim articular cartilage debris around margin of implant. Remove bone particles and lavage thoroughly.

To ensure mechanical interlock of the Taper Post and implant, carefully clean Taper Post taper with provided instruments. All drilling or reaming should be done at slowest speeds possible with vigorous lavage to minimize heat effects to adjacent bone and cartilage tissues.

Accepted practices in post-operative care should be used. The patient is to be instructed and monitored to ensure a reasonable degree of compliance to post-operative instructions and activity restrictions.

Excessive activity, impact, and weight gain have been implicated in the reduction of the benefit and service life of prosthetic devices.

These implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. Their safety in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

#### **Precautions**

These implants are intended to be fitted and installed with the associated instruments. Use of instruments from other systems may result in improper implant selection, fitting, and placement which could result in implant failure or poor clinical outcome. The instruments should be regularly inspected for any signs of wear or damage.

Surgeon or Physician should discuss general risks and potential complications associated with this and any surgical procedure with the patient prior to patient consent.

#### **Possible Adverse Effects**

1. Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions. Particulate wear debris and mild tissue discoloration from metallic components have been noted in other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
2. Infection or allergic reaction.
3. Loosening, migration or loss of fixation of implant.
4. Fretting and crevice corrosion can occur at the interface between the implant components.
5. Fatigue fracture of the implants as a result of bone resorption around the implant components.
6. Wear and damage to the implant articulating surface.
7. Wear and damage to the adjacent and opposed articular cartilage surfaces or soft tissue support structures.
8. Intraoperative or postoperative bone fracture.
9. Postoperative pain or incomplete resolution of preoperative symptoms.

10. Periarticular calcification or ossification, with or without impediment of joint mobility.
11. Incomplete range of motion due to improper selection or positioning of components.
12. Transient nerve palsy.
13. Embolism.

**Sterility**

Implants and single-use disposable instruments are provided STERILE. Metallic implant components are sterilized by exposure to gamma radiation. Non-metallic implant components are sterilized by gas plasma sterilization. The single-use disposable instruments are sterilized by exposure to gamma radiation. Do not resterilize. Do not use components if packaging is opened or damaged. Do not use components if beyond expiration date. Do not reuse implants or single-use disposable instruments. Reuse of these devices can increase the risk of patient infection and can compromise service life and other performance attributes of the device(s).

**Caution**

United States Federal Law restricts this device to sale by or on the order of a physician.

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
1. Use **Drill Guide** to locate the axis normal to the articular surface and central to the defect. Choose the **Drill Guide** to circumscribe the defect. Place **Guide Pin** into a Cannulated Powered Drill and secure at the etch marking on the **Guide Pin**. Advance **Guide Pin** into bone making sure that it is central to the defect. (It is important to verify the **Drill Guide** is seated on the curved surface such that all 4 points of contact are established on the articular surface. A normal axis and correct **Articular Component** diameter are necessary for proper implant fit.)

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
5. Clean taper in **Taper Post** with **Taper Cleaner**. Place **Trial Cap** into **Taper Post** to confirm correct depth of **Taper Post**. The height of the **Trial Cap** must be flush or slightly below the existing articular cartilage surface to avoid the **Articular Component** from being placed proud or above the surface of the defect. Adjust depth if needed using the **Driver** to rotate the **Taper Post** (rotate clockwise to advance and counterclockwise to retract). Remove **Trial Cap**.

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
6. Place **Centering Shaft** into taper of **Taper Post**. Place **Contact Probe** over **Centering Shaft** and rotate around shaft. Use light pressure on the Contact Probe to ensure proper contact with the articular surface. Read **Contact Probe** to obtain offsets at indexing points and mark each of the identified offsets on the appropriate Sizing Card. Select appropriate **Articular Component** using Sizing Card.

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
11. Before placing the **Articular Component** on the **Implant Holder** make sure that sufficient suction is present to hold the device on the distal suction cup. Align the **Articular Component** on the **Implant Holder**. For non-spherical **Articular Components** orient the etch marks on the back of the **Articular Component** with the etch mark on the handle of the **Implant Holder**. Align the **Articular Component** with the appropriate offsets. Insert into taper of **Taper Post**.

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---

### Revision Procedure

In the event that the HemiCAP® requires removal, please follow these steps.

1. Remove any tissue directly on and up to the outer edge of the HemiCAP® device. Choose the appropriate diameter **Revision Cutter** and place into a Cannulated Powered Drill. Advance the **Revision Cutter** counterclockwise through the tissue surrounding the outer edge of the device until all the tissue is removed and the full depth of the outer rim of the HemiCAP® device is completely exposed.
2. Place the appropriate diameter **Revision Driver** over the HemiCAP® device and then impact the **Revision Driver** until it locks around the edge of the HemiCAP® implant. Twist the **Revision Driver** in a counterclockwise rotation until the implant is removed.
3. If necessary, use an osteotome or rongeur to remove the tissue surrounding the head of the Taper Post until sufficient space is available to insert the **Female Hex Driver**. Insert the Driver onto the Taper Post and then rotate counterclockwise to remove the Taper Post from the bone.

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	Consult Instructions for Use
	Caution
	Do Not Reuse
	Sterilized using Radiation
	Sterilized using Vaporized Hydrogen Peroxide
<b>R<sub>x</sub> ONLY</b>	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician
	Do Not Use If Package Is Damaged
	Lot Number
	Expiration Date
	Manufacturer
	Catalogue Number
	Single sterile barrier system with protective packaging inside
	Double sterile barrier system
	Date of Manufacture

Arthrosurface, Inc.  
319 Manley Street,  
West Bridgewater, MA 02379  
tel +1 508 520 3003  
www.arthrosurface.com

**R<sub>x</sub>ONLY**

**STERILE R**

This product is covered by one or more of  
U.S. Patent Nos. 6,520,964; 6,610,067; 6,679,917;  
other patents and other patents pending.  
HemiCAP® is a trademark of Arthrosurface, Inc.  
U.S. © 2024 Arthrosurface, Inc. All rights reserved.  
Printed in the U.S.A.

PN 4001-2002 REV N