Instructions for Use

Total Toe – Proximal Phalanx Articular Component

Description

The Arthrosurface ToeMotion® Total Toe System consists of a metatarsal component and a phalangeal component designed for arthroplasty of the 1st metatarsal head and the base of the proximal phalanx. These two implants replace the metatarsophalangeal joint by complete functional preservation of the joint and maintaining of the sesamoid complex.

The enclosed Phalangeal Fixation Component is to be used with an appropriate HemiCAP® MTP metatarsal head articular implant (sold separately).

Materials

Proximal Phalanx Fixation Component:

Titanium Alloy (Ti-6Al-4V)

Proximal Phalanx Articular Insert:

Ultra High Molecular Weight Polyethylene (UHMWPE)

Sold Separately

Metatarsal Head Articular Component:

Cobalt-Chromium Alloy (Co-Cr-Mo)

CP Ti Plasma Spray Coating applied

Indications

The Arthrosurface Total Toe Arthroplasty System is a two-piece implant that is intended to be used as prosthesis for the metatarsophalangeal joint (MTP). The device is intended for cemented use only. Indications for use include:

- Painful degenerative metatarsophalangeal joint change
- Hallux rigidus stage 3 and 4
- Hallux valgus and hallux rigidus
- Hallux limitus with painful arthrofibrosis
- Revisions after moderate proximal phalanx resection

Patient selection factors to be considered include:

1. Need to obtain pain relief and improve function.
2. Patient age as a potential for early-age revision of total joint arthroplasty.
3. Patient overall well-being, including ability and willingness to follow instructions and comply with activity restrictions.
4. Failure of previous conservative treatment options in correcting deformity and achieving pain relief.

Contraindications

Absolute contraindications include:

1. Significant bone demineralization or inadequate bone stock.
2. Inadequate skin, musculotendinous or neurovascular system status.
3. Inflammatory or rheumatoid arthritis, infection, sepsis, and osteomyelitis.
4. Patients that have a known sensitivity to Cobalt-Chrome alloys or UHMWPE polymers typically used in prosthetic devices.

Relative contraindications include:

1. Uncooperative patient or patient incapable of following preoperative and postoperative instructions.
2. Metabolic disorders which may impair the formation or healing of bone.
3. Infections at remote sites which may spread to the implant site.
4. Rapid joint destruction or bone resorption visible on roentgenogram.
5. Chronic instability or deficient soft tissues and other support structures.
6. Vascular or muscular insufficiency.

Warnings

Improper selection, placement, positioning, alignment, and fixation of the implant components may reduce the service life of the prosthetic components. Inadequate preparation and cleaning of the implant components mating surfaces may result in improper fixation of the device. Improper handling of the implants can produce scratches, nicks or dents that may have adverse clinical effects on mating joint surfaces. Do not modify implants. The surgeon shall be thoroughly familiar with the implants, instruments, and surgical technique prior to performing surgery.

Incorrect sizing of the Proximal Phalanx implant could result in excessive loading of the affected toe joint.

Accepted practices in post-operative care should be used. The patient is to be instructed and monitored to ensure a reasonable degree of compliance to post-operative instructions and activity restrictions.

Excessive activity, impact, and weight gain have been implicated in the reduction of the benefit and service life of prosthetic devices.

These implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. Their safety in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Precautions

These implants are intended to be fitted and installed with the associated instruments. Use of instruments from other systems may result in improper implant selection, fitting and placement, which could result in implant failure or poor clinical outcome. The instruments should be regularly inspected for any signs of wear or damage.

Surgeon or Physician should discuss general risks and potential complications associated with this and any

surgical procedure with the patient prior to patient consent.

Possible Adverse Effects

1. Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions. Particulate wear debris and mild tissue discoloration from metallic components have been noted in other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
2. Infection or allergic reaction.
3. Loosening, migration or loss of fixation of implant.
4. Fretting and crevice corrosion can occur at the interface between the implant components.
5. Fatigue fracture of the implants as a result of bone resorption around the implant components.
6. Wear and damage to the implant articulating surface.
7. Wear and damage to the adjacent and opposed articular cartilage surfaces or soft tissue support structures.
8. Intraoperative or postoperative bone fracture.
9. Postoperative pain or incomplete resolution of preoperative symptoms.
10. Periarticular calcification or ossification, with or without impediment of joint mobility.
11. Incomplete range of motion due to improper selection or positioning of components.
12. Transient nerve palsy.
13. Embolism.

Sterility

Implants and single-use disposable instruments are provided STERILE. Metallic implant components are sterilized by exposure to gamma radiation. Non-metallic implant components are sterilized by gas

plasma sterilization. The single-use disposable delivery tool instrument is sterilized by exposure to EtO Gas. All other single-use disposable instruments are sterilized by exposure to gamma radiation. Do not resterilize. Do not use components if packaging is opened or damaged. Do not use components if beyond expiration date. Do not reuse implants or single-use disposable instruments. Reuse of these devices can increase the risk of patient infection and can compromise service life and other performance attributes of the device(s).

Caution

United States Federal Law restricts this device to sale by or on the order of a physician.

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4. Remove the **1.5mm Guide Pin**, and place the **Fixation Component** into the tapped pilot hole. Using the **Hex Driver**, turn the **Fixation Component** counter clockwise to align threads of **Fixation Component** to the tapped pilot hole. Advance the **Fixation Component** clockwise until the distal surface of the **Fixation Component** is fully seated in the prepared bone socket.

*Note: Before inserting the **Fixation Component**, thoroughly cleanse the pilot hole of any debris and then inject the cement in a retrograde fashion from the end of the hole upwards.*

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	Consult Instructions for Use
	Caution
	Do Not Reuse
	Do Not Re-Sterilize
	Sterilized using Radiation
	Sterilized using Vaporized Hydrogen Peroxide
	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician
	Do Not Use If Package Is Damaged
	Lot Number
	Expiration Date
	Manufacturer
	Catalogue Number
	Single sterile barrier system with protective packaging inside
	Double sterile barrier system
	Date of Manufacture

Arthrosurface, Inc.
319 Manley Street
West Bridgewater, MA 02379
tel +1 508 520 3003
www.arthrosurface.com

R_x ONLY

STERILE | R

This product is covered by one or more of U.S.
Patent Nos. 6,520,964; 6,610,067; 6,679,917;
other patents and other patents pending.
HemiCAP® is a trademark of Arthrosurface, Inc. U.S.
© 2024 Arthrosurface, Inc. All rights reserved.
Printed in the U.S.A.

PN 4001-2005 REV H