UniCAP
Hemiarthroplasty Components
UniCAP® Instructions for Use

Unicompartmental Knee Resurfacing Prosthesis (UniCAP®) Instructions for Use

Description

The Unicompartmental Knee Resurfacing Prosthesis (UniCAP®) incorporates a low-profile femoral articular component that mates to a taper post via a taper interlock. The femoral component articulates against an all-polyethylene tibial component. The UniCAP® implants allow resurfacing of the compartment utilizing the undisturbed compartmental structures and soft-tissues.

Materials

Femoral Components	
Articular Resurfacing Component:	Cobalt-Chromium-Molybdenum alloy (Co-Cr-Mo)
Surface Coating:	Titanium (CP Ti)
Taper Post:	Titanium Alloy (Ti-6Al-4V)
Tibial Components	
Ultra-High Molecular Weight Polyethylene (UHMWPE)	

Indications for Use

Partial replacement of the articulating surfaces of the knee when only one side of the joint is affected due to the compartmental primary degenerative or post-traumatic degenerative disease, previous tibial condyle or plateau fractures, deformity or revision of previous arthroplasty. This device is intended to be used with bone cement.

Patient selection factors to be considered include:

1. Patient's need to obtain pain relief and improve function is significant.
2. Patient's joint stability is good, with intact ACL, and limited mechanical axis deformity.
3. Patient's overall well-being is good, including an ability and willingness to follow instructions and comply with activity restrictions.

Contraindications

Absolute contraindications include:

1. Infection, sepsis, or osteomyelitis.
2. Patients that have a known sensitivity to materials typically used in orthopedic prosthetic devices or bone cements.

Relative contraindications include:

1. Uncooperative patient or patient incapable of following pre-operative and post-operative instructions.
2. Osteoporosis.
3. Metabolic disorders which may impair the formation or healing of bone.
4. Infections at remote sites which may spread to the implant site.
5. Rapid joint destruction or bone resorption visible on roentgenogram.
6. Chronic instability or deficient soft tissues and other support structures.
7. Vascular or muscular insufficiency.
8. Inadequate skin, musculotendinous or neurovascular system status.

Warnings

Improper selection, placement, positioning, alignment, and fixation of the implant components may reduce the service life of the prosthetic components. Inadequate preparation and cleaning of the implant components mating surfaces may result in improper fixation of the device. Improper handling of the implants can produce scratches, nicks or dents that may have adverse clinical effects on mating joint surfaces. Do not modify implants. The surgeon shall be thoroughly familiar with the implants, instruments, and surgical technique prior to performing surgery.

The all-polyethylene tibial component is intended to be used in conjunction with activity restrictions.

When taking readings of articular surfaces, care should be taken to ensure that the distal tip of contact probe is making contact on articular surfaces and free from any soft tissue structures to ensure accuracy. Use light pressure on contact probe to slightly indent articular surface at chosen points, this will ensure selected implant will be recessed just below articular surface at margins of implant.

When placing implant, carefully trim articular cartilage debris around margin of implant. Remove bone particles and lavage thoroughly. To ensure mechanical interlock of the taper post and implant, carefully clean taper post with provided instruments. All drilling or reaming should be done with vigorous lavage to minimize heat effects to adjacent bone and cartilage tissues.

Ensure that care is taken to obtain complete and uniform bone cement coverage at implant site. Unsupported components or unevenly supported components may result in implant failure.

Accepted practices in post-operative care should be used. The patient is to be instructed and monitored to ensure a reasonable degree of compliance to post operative instructions and activity restrictions. Excessive activity, impact, and weight gain have been implicated in the reduction of the benefit and service life of prosthetic devices.

Precautions

The Unicompartamental Knee Resurfacing Prosthesis (UniCAP®) implant is intended to be fitted and installed with the Unicompartamental Knee Resurfacing Prosthesis (UniCAP®) instrument set. Use of instruments from other systems may result in improper implant selection, fitting, and placement which could result in implant failure or poor clinical outcome. The instrument set should be regularly inspected for any signs of wear or damage. Do not reuse implants. Reuse of single use devices can increase the risk of patient infection and can compromise service life and other performance attributes of the device.

Possible Adverse Effects

1. Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions. Particulate wear debris and mild tissue discoloration from metallic components have been noted in other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
2. Infection or allergic reaction.
3. Loosening, migration or loss of fixation of implant.
4. Fretting and crevice corrosion can occur at the interface between the implant components.
5. Fatigue fracture of the implants as a result of bone resorption around the implant components.
6. Wear and damage to the implant articulating surface.
7. Wear and damage to the adjacent and opposed articular cartilage surfaces or soft tissue support structures.
8. Intraoperative or postoperative bone fracture.
9. Postoperative pain or incomplete resolution of preoperative symptoms.
10. Periparticular calcification or ossification, with or without impediment of joint mobility.
11. Incomplete range of motion due to improper selection or positioning of components.
12. Transient peroneal palsy.

Sterility

Metallic prosthetic components are sterilized by exposure to gamma irradiation. Non-metallic prosthetic components are sterilized by gas plasma sterilization. Do not resterilize any components. Do not use components if packaging is opened or damaged. Do not use components if beyond expiration date.

Caution: United States Federal law restricts this device to sale by or on the order of a physician.

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
1. The tibial component may be placed under arthroscopic visualization to ease preparation of the tibial defect. With the knee at 90 degrees of flexion, and working through an anteromedial portal, place a **Tibial Template** central to the damaged area of the tibial plateau surface. The underside of the **Tibial Template** should be in contact with the surface of the tibia. Select the **Tibial Template** that best matches the A/P and M/L surface curvatures of the tibia. Attach the **Tibial Template** to the **Tibial Drill Guide** and **Bullet**. Place the guide so the arm of the **Tibial Template** is parallel to the tibial plateau.

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2. Drill the **2.0mm Drill Tipped Pin** through the central axis of the **Tibial Drill Guide** until it reaches the center of the **Tibial Template**. Use a small closed curette to “catch” the tip of the **Pin** to prevent drilling into the femur. Care must be taken so that excessive torque is not applied to the **Drill Guide** which may cause the pin to miss the target. Confirm that there is a minimum of 5mm of bone from the edge of the **Template** to front of the tibia to avoid breaking through the anterior tibia during reaming. Bring the knee into extension and probe the underside of the meniscus to ensure complete visualization and proper placement of the **Template**. Remove the **Tibial Drill Guide** and **Bullet**.

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4. Advance the **Introducer** into the prepared tibial tunnel. The proximal tip of the **Introducer** should be flush with the tibial plateau. Begin to advance the threaded **Blade Stop** over the **Introducer** until it begins to screw into the bone. Remove **Introducer** and continue to advance **Blade Stop** until it is 2/3rds into the tunnel.
 - a. Remove **Driver** handle and reinsert **Introducer** and **Driver** and continue to advance as one unit. Stop when the tip of the **Introducer** is flush with the tibial plateau. Confirm that the laser mark on the **Introducer** is in-line with the laser mark in the slotted window of the **Blade Stop Driver**.
 - b. **Blade Stop** is at the correct depth when the tip of the **Introducer** is flush with the tibial plateau and the laser mark lines on the **Driver** and **Introducer** are aligned in the slotted window. Remove the **Blade Stop Driver** and **Introducer**.

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
5. Place the **Cutting Blade** into the **Blade Holder** with the long slot facing posteriorly to the joint. Introduce the **Cutting Blade** into the portal. Advance the **Blade Drive Shaft** into the tibial tunnel until it is visible in the joint. Push the tip of the **Drive Shaft** through the center of the **Cutting Blade**.

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6. Connect the **Cutting Blade** to the **Blade Drive Shaft** by turning the **Drive Shaft** 90 degrees and pull distally to engage the blade.
 - a. To lock the Tibial Cutting system, push the sheath in an upward motion and rotate 90 degrees so the Lock Indicator on **Blade Drive Shaft** is positioned over the Dowel Pin and Laser Mark line. Release to lock **Cutting Blade** into position. Attach the Powered Drill to the laser mark indicated on the distal end of the **Drive Shaft**.

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7. Using a Powered Drill, begin rotation counterclockwise to normalize the blade to the tibial plateau. This will help ensure even cutting engagement of the **Cutting Blade** into the plateau. Care should be taken to avoid the meniscus. Begin to prepare the

inlay implant socket using a clockwise blade rotation. Drilling is complete when the **Cutting Blade** comes in contact and stops on the proximal end of the **Blade Stop**.

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
8. With **Cutting Blade** remaining in place, insert the appropriate sized **Sizing Trial** based on the mapping determined with the **Tibial Template**. Confirm fit at anterior/-posterior and medial/lateral margins of the

Sizing Trial. If the **Trial** is proud at the margins, use the **Blade Stop Wrench** to rotate the **Blade Stop** clockwise. Each 90 degree turn of the **Blade Stop** with **Blade Stop Wrench** lowers the **Blade Stop** 1mm. Reattach the Powered Drill and re-ream to the new depth. This will lower the implant depth in the tibial socket.

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9. With the **Sizing Trial** set at the appropriate height, begin removing the instrumentation. Raise the **Drive Shaft** so the **Cutting Blade** can be grasped for removal. To unlock and remove the **Cutting Blade** push the sheath in an upward motion and rotate 90 degrees counter clockwise so the Unlock Indicator on **Blade Drive Shaft** is positioned over the Dowel Pin. Release to unlock **Cutting Blade**. Push **Drive Shaft** upwards to free it from the **Cutting Blade** and turn the **Drive Shaft** 90 degrees. Pull distally on the **Drive Shaft** and remove from the tibial tunnel. Use a grasper to remove the **Cutting Blade**.

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
12. Implant **Tibial Component** with cement using the **Arthrosurface® Cement Ejector**. When discharging cement, allow the back pressure from the cement extraction to lift the implant up 2mm and then continue to back fill the tibial tunnel with cement. Utilize the **Tibial Template** through the portal to apply downward pressure onto the **Tibial Component** to seat it in its final position. This will allow for optimal cement integration.

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
Cement Ejector Assembly

- Mix low viscosity cement according to manufacturers' directions.
- Place cement into the **Delivery Syringe**. Remove **Funnel** when **Delivery Syringe** is full.
- Insert **Plunger** into **Delivery Syringe**. Insert into **Ejector Handle**.
- Attach **Drive Rod** onto **Ejector Handle**.
- Place **Threaded Sheath** into tibial tunnel to prevent cement extrusion.
- A powered drill with Jacobs chuck is used to advance the **Drive Rod**.
- Deliver cement to undersurface of implant through tibial tunnel and draw **Cement Ejector** retrograde as tibial tunnel is filled.

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
2. Drill the **2.0mm Threaded Pin** through the central axis of the **Femoral Drill Guide** into bone to the laser mark line on the **2.0mm Threaded Pin**. Remove the **Femoral Drill Guide**.

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
3. Drive the **Femoral Centering Shaft** over the **Threaded Pin** so the laser mark line is at the height of the original articular surface. Good practice is to leave the etch mark proud initially and the advance by hand with the **Driver**.

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
4. Place the 40mm **Contact Probe** over the **Femoral Centering Shaft**. Read the **Contact Probe** to take the superior and inferior offsets and mark them onto the appropriate sizing card.

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5. Repeat using the 20mm **Contact Probe** to obtain the offsets medially and laterally and mark the sizing card.

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
6. Select the appropriately sized **Central Femoral Reamer** based on the average medial to lateral mapped offset. This will either

be a 2mm or 3mm reamer. Prepare central femoral cut by advancing the **Central Femoral Reamer** over the **Centering Shaft** until it contacts the stop. Remove **Centering Shaft**.

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
7. Select the appropriately sized **Guide Block** based on the average anterior/posterior offset from the sizing card (6mm to 10mm by 1mm increments) and attach it to the **Femoral Drill Guide**. Realign **Femoral Drill Guide** on the distal femur. Maintain four points of contact to ensure accurate placement of guide pins.
 - a. Insert the **Pin Sleeves** into the slots located superior and inferior on the **Femoral Drill Guide**. Beginning with the superior **Pin Sleeve**, drill the **Short Threaded Pins** through the **Pin Sleeves** into bone to the laser mark line on the **Short Threaded Pin**. Remove the **Pin Sleeves** and remove the **Femoral Drill Guide**.

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
10. Select the appropriately sized **Femoral Sizing Trial**, attach to **Sizing Trail Handle** and place into position. Leave superior and inferior **Short Threaded Pins** in place. Confirm fit at anterior/posterior margins and medial/lateral margins of the **Femoral Sizing Trial**.

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
11. Before preparing the pilot hole for the **Taper Post**, be sure **Femoral Trial** is seated so the edges are flush or slightly recessed to the cartilage. Advance the **Femoral Pilot Drill** through the **Sizing Trial Handle** until the laser mark is flush with the end of the handle and leave it in position. Replace the **2.0mm Threaded Pin**.

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
14. Insert the **Taper Post** into the **Sizing Trial Handle**. Attach the **Handle** and **Taper Post** assembly onto the **Femoral Sizing Trial**. Insert the **Hex Driver** into the **Handle** and advance the

Taper Post into the bone. Stop advancing the **Hex Driver** when the raised stop on the **Driver Shaft** contacts the top of the **Sizing Trial Handle** and the **Femoral Sizing Trial** is flush with the surrounding cartilage. Remove the **Sizing Trial**.

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15. Use the **Final Placement Gauge** in combination with the **Femoral Trial** to confirm the proper depth placement of the **Taper Post**. With the **Femoral Trial** in position, and with the **Final Placement Gauge** coupled to the **Taper Post**. A very small separation (less than .5mm) should be seen between the components. If no gap is visible, the **Taper Post** is set too deep and should be raised. If a large gap is visible, the **Taper Post** is set too shallow and should be lowered. Be sure not to impact the **Final Placement Gauge** onto the **Taper Post**.

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
	Consult Instructions for Use
	Caution
	Do Not Reuse
	Do Not Re-Sterilize
	Sterilized using Radiation
	Sterilized using Vaporized Hydrogen Peroxide
R_x ONLY	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician
	Do Not Use If Package Is Damaged
	Lot Number
	Expiration Date
	Manufacturer
	Catalogue Number
	Single sterile barrier system with protective packaging inside
	Double sterile barrier system
	Date of Manufacture

Manufacturer

Arthrosurface, Inc.
319 Manley St, W Bridgewater, MA 02379
tel +1 508 520 3003 • fax +1 508 528 3785
www.arthrosurface.com

STERILE	H ₂ O ₂
---------	-------------------------------

STERILE	R
---------	---

R_xONLY

This product is covered by one or more of U.S. Patent Nos.
6,520,964; 6,610,067; 6,679,917; other patents and other patents
pending. HemiCAP® is a trademark of Arthrosurface, Inc. U.S.
© 2019 Arthrosurface, Inc. All rights reserved. Printed in the U.S.A.

PN K001-2001 REV L