



# ToeMate

## Hammertoe Correction System

### A Strong Step in the Right Direction

ToeMate® is an easy to implant bone screw system intended for the correction of hammertoe deformity. It is provided in a complete packaged and sterile kit.

- Proven, strong & stable fixation
- Quick & accurate joint alignment
- Disposable instrument kit



#### Available in 3 Sizes:

- Large - Ø 4.7mm
- Small - Ø 3.8mm
- X-Small - Ø 3.2mm

*\* Implants not shown to scale*



Implant components completely embedded in the phalangeal bones (i.e. no implant exposure post-op) avoids complications associated with pin tract infection.



All instruments hand driven. No power tools or connection required.

*\*Except for x-small*



All implants and instruments provided in a complete packaged & sterile kit. This eliminates missing instruments, incomplete trays and minimizes the risk of infection.

### Ordering Information

#### Implant Package (LARGE - Ø 4.7mm)

9H00-T010	Screw, Proximal (Purple), L: 11mm Screw, Distal (Blue), L: 6.5mm Taper Lock Pin (0° & 10°) Delivery Tool (Bone White)
-----------	--

#### Implant Package (SMALL - Ø 3.8mm)

9H00-T030	Screw, Proximal (Purple), L: 9mm Screw, Distal (Gold), L: 7mm Taper Lock Pin (0° & 10°) Delivery Tool (Blue)
-----------	---

#### Implant Package (X-SMALL - Ø 3.2mm)

9H00-T050	Screw, Proximal (Purple), L: 9mm Screw, Distal (Gold), L: 7mm Taper Lock Pin (0° & 10°) Delivery Tool (Teal)
-----------	---

#### Instrument Package (LARGE - Ø 4.7mm)

9H00-T020	Drill Screw Driver (Bone White) Thimble
-----------	---

#### Instrument Package (SMALL - Ø 3.8mm)

9H00-T040	Drill Screw Driver (Blue) Thimble
-----------	---

#### Instrument Package (X-SMALL - Ø 3.2mm)

9H00-T060	1.5mm Guide Pin Screw Driver (Teal) Thimble
-----------	---

**Arthrosurface, Inc.**

319 Manley Street, West Bridgewater, MA 02379  
508-520-3003 | [movefreely@arthrosurface.com](mailto:movefreely@arthrosurface.com)

[www.arthrosurface.com](http://www.arthrosurface.com) | Arthrosurface. Restore. Recover. Renew. | Move Freely | Preserve Motion

This document and information is intended for markets where regulatory approval has been granted.

Arthrosurface, ToeMate, Move Freely, Preserve Motion, and Restore. Recover. Renew. are trademarks and/or registered trademarks of Arthrosurface, Inc. and its affiliates in certain jurisdictions.

System designed and manufactured in the USA | Printed in the USA

©2025 Arthrosurface, Inc. All rights reserved.

PN 4001-4008 REV E

