

WristMotion[®]

Total Wrist Arthroplasty System

Designed to preserve motion and maximize carpal stability, the WristMotion Total Wrist Arthroplasty (TWA) System provides an advanced solution for wrist arthritis

Secure Carpal Fixation¹


- Rock-solid fixation with the novel center threaded taper post
- Preserves anatomic center of rotation on the capitate by adjusting height on the radial side
- Proven Morse Taper interlock for enhanced rotational stability
- Reduces torque on carpal components by adjusting the tension through radial offset heights²

Improved Range of Motion and Rotational Freedom¹

- Greater range of extension due to the increased surface area of the dorsal flange
- Enables 20 degrees of axial rotational from the radial poly
- Maintains Dart Thrower's Motion

Precise Anatomic Fit¹

- The guided instrumentation enables precise and consistent implant placement



The Next Generation of Wrist
Arthroplasty

WristMotion *Total Wrist Arthroplasty System*

Personalized implants for varying patient anatomies


- Radial and Carpal Implants available in 3 sizes (S/M/L)
- Radial Stems with Ti-Plasma Spray coating available in left and right options
- Radial Polys available in 3 height offsets (+0mm, +1.5mm, +3mm)
- Bone Screws available in lengths 10mm to 45mm



Individualized Care Through a Comprehensive Portfolio

Our full wrist arthroplasty portfolio provides flexibility to select the right product and sizes for the patient

WristMotion® Hemiarthroplasty Implant System



WristMotion® Total Wrist Arthroplasty System



*X-ray includes 3.0mm radial poly

More than 2,600 WristMotion® Hemiarthroplasty implants performed^{1,3}

For complete product information, please visit: www.arthrosurface.com

1. Preclinical data on file. Results may not correlate to clinical performance
2. Reference: Akhbari, B. et al. Proximal-distal shift of the center of rotation in a total wrist arthroplasty is more than twice of the healthy wrist. Journal of Orthopedic Research, Vol. 28, Issue 7, p. 1575-1586
3. As of July 31, 2021. Data on file.

Arthrosurface. Restore. Recover. Renew. | Move Freely | Preserve Motion

319 Manley Street, West Bridgewater, MA 02379 | 508-520-3003 | movefreely@arthrosurface.com

©2025 Arthrosurface, Inc. | PN 3001-4016 REV B

Arthrosurface and WristMotion are registered trademarks of Arthrosurface


arthrosurface®